



This sponge

By eliminating the disadvantages of earlier methods, the Triosorb Sponge has achieved a real breakthrough in thyroid testing. **It is an in vitro test unmatched in accuracy, speed and convenience.**

Accuracy: Because factors such as red blood cells and exogenous iodine have been eliminated from consideration in the Triosorb Test, it is unmatched in accuracy.

Speed: With only 3 washes and no need for double pipettings, shakers, or incubators, the Triosorb Test can be

revolutionized thyroid testing!

more rapidly performed than any other T-3 test.

Convenience: Triosorb is in a disposable kit ready for immediate use at room temperature, making it the simplest and most convenient thyroid function test to perform.

McAdams* reported that "The resin sponge (Triosorb) technique is superior to the erythrocyte method for performing the I^{131} T3 test in terms of simplicity, convenience and elimination of errors characteristic of the erythrocyte procedure."

Triosorb is available to all doctors, hospitals and clinical laboratories—AEC licensing is not required. Because Triosorb will enable far more screenings to be performed, this procedure may soon become as standard as today's blood counts and urinalyses.

801202




*McAdams, G. B. and Reinfrank, R. F., *Jrnl. Nuclear Med.*, 5:112, Feb., 1964.

TRIOSORB[®]
T-3 DIAGNOSTIC KIT
ABBOTT LABORATORIES NORTH CHICAGO, ILL.

RADIOISOTOPE SCANNER

MODEL 54-FD

DUAL, OPPOSED, 5-INCH CRYSTALS




The demonstrable advantages of a dual 5-inch crystal scanner should be investigated by all those with a high clinical load who desire high resolution, rapid scans of both large and small organs or of the whole body.

The two scanning heads, exactly opposite each other, have separate, and complete electronics and print-out so that the data collected by each crystal may be used separately, in coincidence, or additively.

Mechanical and electronic specifications are the same as for our other large-crystal radioisotope scanners Models 54F and 54H:

Scanning speeds continuously variable to 200 inches per minute (500 cm/min.); adequate shielding even for high energy gamma emitters (up to 3 inches lead and 1 inch steel); high resolution crystals (9 per cent or better); accurate, reproducible scanning speeds and line spacing; no scalloping at any speeds; low background crystals (2 inch thick pure NaI light pipe); Gamma-graphic (patent pending) or slit mask photoscans; unequivocal one year warranty anywhere in USA or Canada.

This unparalleled radioisotope scanner is priced at \$28,750 with delivery in 90 days guaranteed.




OHIO-NUCLEAR, INC.

1725 FALL AVENUE

CLEVELAND, OHIO

216 - 621-8477




SINGLE SOURCE

For the New Short
Half-Life or Low Energy
Isotopes


More and more users are looking to Nuclear Consultants Corporation for the latest developments in radiopharmaceuticals. Whatever your need, whenever it arises, you can depend on NCC for demonstrated leadership and service.

Chlormerodrin-Hg¹⁹⁷


Chlormerodrin-Hg²⁰³

 Proven brain and renal scanning agents. Hg¹⁹⁷ available at lower cost through NCC's "Controlled Decay" program.

1-Mercuri-2-Hydroxypropane-Hg¹⁹⁷

 For improved spleen scanning. Much simpler and more convenient than heat treating tagged red cells. MHP-Hg²⁰³ also available for experimental use.


Iodine-I125

 I-125 is a longer-life isotope than I-131, and offers significantly lower radiation dosage to the patient. Available as tagged albumin and other iodinated compounds.


Technetium-99m

For Scanning Brain...Thyroid...Liver...Lung


TechnKow[®] generator

 Supplies ample material as sodium pertechnetate for daily use. Economical for as few as seven or eight brain scans per week. New #1810 Shielded Dispenser available — simple semi-automated method of obtaining sterile solution.

"Instant Technetium"

 Sodium pertechnetate in sterile, pyrogen-free solution permits immediate use without the necessity for "milking," sterilization, or calibration.

Technetium-Sulfide Colloid

 In sterile, pyrogen-free form, reported as providing excellent liver scans. Also usable for lung scanning following the technique of Quinn. 6517



**NUCLEAR CONSULTANTS
CORPORATION**


Box 6172, Lambert Field • St. Louis, Mo. 63145
314 PErshing 8-8927

CHICAGO • CLEVELAND • HOUSTON • LOS ANGELES • NEW YORK • SAN FRANCISCO • WASHINGTON, D.C.

FOR BIOLOGICAL RESEARCH

... CLINICAL DIAGNOSIS

... SPECIAL STUDIES



ARMAC® provides rapid, accurate measurement of radioactivity

The Armac Scintillation Detector is ideally suited for *in vivo* research using gamma-emitting tracers. Its 4-pi counting geometry and high counting efficiency require the administration of only small doses of radioactivity, yet assure significant measurements in short counting periods.

- In biological research Armac can be used to count radioactivity in small laboratory animals in the study of retention and excretion of various gamma-emitting isotopes. Because only small doses of radioactivity are required, sacrificing is unnecessary and the animal acts as its own control for repeated experimentation.
- In clinical research Armac is useful in nutritional

studies and other applications where changes in the amount of circulating radioactivity in the blood can be directly related to other body functions. These changes can be measured accurately and rapidly using the blood-rich portion of the patient's forearm positioned in the Armac's counting chamber.

- In special studies Armac can be used to quickly and accurately determine the presence of low levels of radioactivity in bulk samples of meat, milk, water, blood, wastes, soil, or tissue with little or no sample preparation. The photograph at right, above, shows an Armac Detector counting radioactivity in a live fish in flowing stream water. (Courtesy Oak Ridge National Laboratory.)

A complete counting system includes the Model 446 Armac Scintillation Detector and any of several Packard Spectrometers, control units and ratemeters. There are more than ten Armac Detector Systems available to meet your specific counting requirements. Choose visual or recorded data output; digital or analog data presentation. Any system can be readily modified or expanded to meet changing research requirements. Special systems incorporating automatic sample changers are also available.

For detailed information contact your Packard Sales Engineer, or write for Bulletin 1024.



PACKARD INSTRUMENT COMPANY, INC.
2200 WARRENVILLE ROAD • DOWNERS GROVE, ILLINOIS 60515
TELEPHONE: 312/969-6000

ANNOUNCING NEW REPRINTS

Journal of Nuclear Medicine

Vols. 1-4. Chicago 1960-1963

Available 1966

Cloth bound set\$90.00
Paper bound set 80.00
Per volume, paper bound 20.00

Atompraxis

International Monthly for Applied Atomic Energy in Industry, Agriculture, the Natural Sciences, and Medicine, with Special Consideration of Medical Radiation Biology and Radiological Safety.

Now Available

Vol. 1. Karlsruhe 1955

Paper bound volume\$10.00

Vols. 2-3, 1956-1957

Per volume, paper bound 20.00

Nuclear Science Abstracts

Now Available

NUCLEAR SCIENCE ABSTRACTS has gained international recognition as the only abstract journal devoted exclusively to announcing, abstracting, and indexing the literature of nuclear science and technology.

Vols. 1-8. Washington 1948-1954

Cloth bound set\$230.00

Paper bound set 200.00

Vols. 1-3, 1948-1949

Per volume, paper bound 20.00

Vols. 4-8, 1950-1954

Per volume, paper bound 30.00

Vol. 11, 1957

Paper bound in 2 parts (2,204 pp.) 60.00



JOHNSON REPRINT CORPORATION

111 Fifth Avenue, New York, New York 10003

JOHNSON REPRINT COMPANY LTD.

Berkeley Square House, London W.1, England

when investing in a scintillation scanner

what should you look for?

among other things . . .

dependable local service



No matter how good your scanner may be, it's no good to you at all if it's out of commission. Your busy schedule simply can't tolerate the disruptions growing out of equipment "downtime".


Across the United States, there are over 400 Picker-employed servicemen working out of 112 Service Centers and backed up by over 40 service specialists. This service force is substantially larger than that of all other nuclear instrument manufacturers *put together*.

PICKER
nuclear

Magna Scanner

the *versatile scanner* / the *proven scanner*

PICKER NUCLEAR
DIVISION / PICKER X-RAY CORPORATION
WHITE PLAINS, NEW YORK



Photostatic of lungs of female patient, aged 50, showing pulmonary emboli.

ALBUMOTOPE[®]-LS

Squibb Aggregated Radio-Iodinated (¹³¹I) Albumin (Human)

FOR EARLIEST EVIDENCE OF PULMONARY EMBOLISM...

A unique diagnostic feature of the lung scan technique that employs Albumotope-LS as the tracer is its capacity to display ischemic areas produced by small pulmonary emboli before evidence of infarction or ischemia is demonstrated in standard chest x-rays.

This method has been found to be a valuable adjunct to selective pulmonary arteriography with the advantages of being readily available, reliable, simple, direct and free from morbidity.

Albumotope-LS meets all the requirements for the ideal radiopharmaceutical for the scanning of pulmonary infarcts as described by Quinn et al.*

- Exhibits a marked difference in concentration between the infarcted and the uninvolved lung tissue.
- Because concentration in the lungs is so selective, it presents a low radiation background when other organs are compared to the lungs.
- Clears rapidly from the blood stream if the infarct concentration of radioactivity is lower than that in the blood.
- Has a suitable gamma energy for scanning and delivers a low radiation dosage to the patient.
- It is non-toxic.

SQUIBB



Squibb Quality—the Priceless Ingredient

Dosage and Scanning Procedure: Recommended scan doses of 150 to 300 microcuries of aggregated radioalbumin depending on the instrumentation available and the technics employed. Scanning immediately follows administration of intravenous injection. Patient may be placed in a prone or supine position.

Side effects and Precautions: There have been no reported cardiovascular or other untoward effects attributable to Albumotope-LS. The product appears to possess no antigenic properties. One patient with a known history of angioneurotic edema, who had been given Lugol's solution in conjunction with aggregated radioalbumin similar to Albumotope-LS, developed urticaria. Extensive clinical use of Albumotope-LS has not borne out the hypothetical


possibility that particles of large size might induce deleterious cardiovascular or cerebrovascular effects.

Radioisotopes should not be used in pregnant women, nursing mothers, or in patients under 18 years of age unless indications are very exceptional.

Available: As a sterile, pyrogen-free, aqueous suspension. Each cc. contains approximately 1 mg. aggregated human serum albumin labeled with 800-1500 microcuries of Iodine-131 at time of manufacture. Also contains 0.9% benzyl alcohol as a preservative.

*Quinn, J. L., III, Whitley, J. E., Hudspeth, A. S. and Watts F. C.: An approach to the scanning of pulmonary infarcts, *J. Nuclear Med.* 5:1, 1964.

Tracerlab Film Badge Service



where people count!

A dosimetry badge is a personal thing; the health, the very life of the wearer may depend on it. That's why Tracerlab people give so much personal care and attention to each and every film badge that passes through their hands, before and after computer processing and quick return to you. And more dosimetry badges pass through their skilled hands than any other monitoring service. Yours should be among them.

You can also count on Tracerlab people for the safest, most reliable radioactive sources in the world — a wide selection of stock or custom types for standardization and analysis. For dependable products and services in the life sciences and health physics, including a broad range of radioanalytical services, come to Tracerlab — where people count!



TRACERLAB

A Division of Laboratory For Electronics, Inc.
WALTHAM, MASSACHUSETTS 02154

Film Badge Service • Health Physics • Bioassays • Sources • Nuclear
Instrumentation • Radiochemicals • Radioactive Waste Disposal
• Radiation Monitoring Instrumentation • Isotope Applications

Advertising Index

Journal of Nuclear Medicine

December, 1965

Abbott Laboratories	
North Chicago, Illinois	IFC, i
Havas International Service	
Amsterdam-Z, The Netherlands	xii
Johnson Reprint Corporation	
New York 3, New York	vii
Nuclear-Chicago	
Des Plaines, Illinois	BC
Nuclear Consultants Corporation	
St. Louis, Missouri	v, IBC
Ohio-Nuclear	
Cleveland, Ohio	iii
Packard Instrument	
Chicago, Illinois	vi
Picker X-Ray Corporation	
White Plains, New York	viii
Squibb, E. R. & Sons	
New York, New York	ix
Tracerlab	
Waltham, Massachusetts	x, xiv


duphar nuclear pharmaceuticals


Why Holland?

Why should you look for isotopes in Holland? Well, because in Holland you will find a company, specialized in isotope production.

The name is Philips-Duphar. Twelve years of isotope production guarantee the kind of experience you can rely on: Nuclear Pharmaceuticals for medicine... C14/H3 compounds and radiochemicals for research... radiographic sources for industry... Representatives in 70 countries and airway despatch make Philips-Duphar your next-door supplier. Write for complete documentation and price lists.

**N.V. PHILIPS-DUPHAR
AMSTERDAM HOLLAND**
Apollolaan 151
Telegrams:
Vitamine-Amsterdam
Telex: 14232





wishing well?

get your wish
for a **WELL**
counting system
at much less cost

Feature by feature, the Tracerlab well counting system matches or exceeds the specifications of other systems costing up to twice as much. See for yourself: 525 lbs. of lead shielding • modular construction for mobility • a full 3" of lead shielding at the side • 5" of top shielding • hinged-lid top accommodates 5 and 15 cc test tube with lid closed — even 4 1/2" diameter beaker • well shield complete as described: \$1295 — the only one at this price. Extra heavy duty units...\$2,400.

Choose from three different scalars, each one the top-performing unit in its class. Select from two scintillation probes with a full two or three inch diameter crystal, each one with a choice of two different diameter wells.

Well counting systems are becoming indispensable in the modern medical diagnostic lab. It makes sense to settle for the highest performance at the lowest price, from the first name in nuclear counting:



TRACERLAB

A Division of Laboratory For Electronics, Inc.

WALTHAM, MASSACHUSETTS 02154

Richmond, California • Houston, Texas • Malines, Belgium • Sales Offices in Principal Cities
Film Badge Service • Health Physics • Bioassays • Sources • Nuclear Instrumentation • Radiochemicals
Radioactive Waste Disposal • Radiation Monitoring Instrumentation • Isotope Applications

Demonstrated Leadership from the specialist

The TECHNEKOW® Shielded Dispenser

A NEW COMPACT "HOT LAB"

by NCC

A Shielded Self-Contained System
for Producing Sterile Technetium
from the TechneKow Source Generator



Completely lead shielded, the new E-1810 Dispenser provides a semi-automatic system for producing a sterile solution of sodium pertechnetate with minimum radiation exposure. It eliminates the need for awkward shielding and handling devices, and enables the smaller laboratory to use technetium on a daily basis.

Insert generator into shielded upper chamber of Dispenser, where it may be stored for daily use. Each day the internal disposable processing parts are replaced with a fresh sterile set.

Operate The TechneKow Shielded Dispenser In 3 Easy Steps

1 □ Pour eluting solution into generator and replace shielded cap.

2 □ When solution has passed through the TechneKow, open dispenser door, reach in and lift up processing unit to activate operation. The solution automatically passes through a sterilizing filter and is injected into a sterile, pyrogen-free bottle.

3 □ Remove bottle with its shielded jacket. Solution is ready for calibration.

The NCC TechneKow Shielded Dispenser is of heavy welded construction with polished chrome plating. A two-inch thick lead shield surrounds the TechneKow source generator. Lead shielding on the walls and door of the lower processing chamber keeps radiation at a minimum. Disposable processing parts are available in kit form.

Write today for further information.




**NUCLEAR CONSULTANTS
CORPORATION**

Box 6172, Lambert Field • St. Louis, Mo. 63145 • 314 Pershing 9-8927


CHICAGO
CLEVELAND
HOUSTON
LOS ANGELES
NEW YORK
SAN FRANCISCO
WASHINGTON, D.C.

PhoGamma derives from the work of H.O. Anger.

This is PHO/DOT™



This is PHO/GAMMA™



NUC-D-4-245

WHY BOTH?

THE SCANNER VS. THE CAMERA. Which of the two do you need—for diagnosis—for clinical research? Which is better suited to your specific application and work load?

WHY PHO/DOT? Because it's the most advanced isotope scanner available. With human-engineered design for simple, foolproof operation. A one-to-one dot or photo record unmatched for readability and resolution. Growing acceptance: one of every two scanners sold today is a Pho/Dot.

WHY PHO/GAMMA? Because it's the only field-tested, performance-proved scintillation camera. Depicts dynamic body processes through rapid-sequence, stop-motion "isotope movies." Produces and records gamma images three to ten times faster than a scanner (even Pho/Dot). Especially economical for hospitals with heavy scanning loads.

THE FACTS ON BOTH? From your Nuclear-Chicago sales engineer—and from your colleagues who are already taking advantage of these remarkable instruments. Or write us.



NUCLEAR-CHICAGO

A DIVISION OF NUCLEAR-CHICAGO CORPORATION

313 Howard, Des Plaines, Ill. 60018 U.S.A.

In Europe: Donker Curtiusstraat 7
Amsterdam W, The Netherlands

Scientists and engineers interested in challenging career opportunities are invited to contact our personnel director.